

Q.Light® 70 NT HOME CARE UNIT Light Therapy for medical and cosmetic application at home


The modern Q.Light® Light Therapy can be used for various indications as standalone or as a ideal complimentary treatment to other methods. It is also applied for **prevention** and **support of healing processes**.

The treatment with Q.Light® **HOME CARE UNIT** supports the different phases of the body's own regeneration and healing potential. It acts inflammation retarding and reduces pain, activates the Immune defense, supports the development of new blood vessels and optimizes microcirculation, helps lymphatic flow and triggers metabolic processes. The Q.Light® Light Therapy regulates the hormonal system and optimizes die neuroimmunological parameters.

The light spectrum of the Q.Light® **HOME CARE UNIT** is clearly defined, without UV, with continues power density and a spectrum of 385 to 780 nm. The emitted light is polarized and has a power density, at max. treatment distance, of 40 mW/cm². The light penetrates deep into skin and tissue. Die dose applied can be individually adjusted.


The Q.Light® Light Therapy can be optimally combined with other standard forms of therapy to enhances notably regeneration capacity.

Side effects have not been reported.



Q.Light® **HOME CARE UNIT**
on **STAND-Home**

The Q.Light® **HOME CARE UNIT** is specially designed for medical and cosmetic application at home. It complies with the Medical Product Regulatory based on the guideline 93/42 EWG.



Manufacturer: Q.Products AG; Schlossgut Wigen; CH - 9404 Rorschacherberg

Fon: +41 (0) 71 858 20 60, Fax: +41 (0) 71 858 20 61; Email: contact@QLight.info - www.QLight.info

Q.Light® 70 NT HOME CARE UNIT

Details, Accessories and Technical Data

Q.Light® HOME CARE UNIT offers, through the emitted full spectrum, a wide range of applications. The units and is easy to handle for anyone .

The special, patented light source produces a pure white light following applications:


Medicine:

- Ø Skin irritations
- Ø Wounds
- Ø Burns
- Ø Inflammations
- Ø Depression
- Ø Sleep disorders
- Ø Energy deficiency
- Ø Strengthening the immune system

Cosmetic:

- Ø Skin reddening
- Ø Wrinkles
- Ø Coupe rose
- Ø Improvement of complexion
- Ø Acne
- Ø Cellulite
- Ø In-depth activity of cosmetic products



Technical Data Q.Light® 70 NT HOME CARE UNIT	
Current	230 V, 50 – 60 Hz
Power supply	50 VA max.
Safety Class	II, Type B
Ø Power density	40 mW/cm ² , 2,4 Joule/cm ² (min.) at max. treatment distance
Spectrum	385 – 780 nm 
Polarization	98%
Treatment area	Ø 25 cm fixed
Emission	Visible light between 385 – 780 nm
Infrared Emission	no IR
UV Emission	no UV
Weight	1'120g
Size	260 x 158 x 173 mm (L x B x H)
Guarantee	24 Month after delivery
Certification	ISO 9002 & EN 46002
CE Mark	1275
Stands	Stand Home, Stand Standard, Stand-Wagen PRO

Manufacturer: Q.Products AG; Schlossgut Wigen; CH - 9404 Rorschacherberg

Fon: +41 (0) 71 858 20 60, Fax: +41 (0) 71 858 20 61; Email: contact@QLight.info - www.QLight.info